

PharmaFocus™ Tubing



HelixMark PharmaFocus™ Tubing

The Right Testing.

The platinum-cured silicone elastomer used for our PharmaFocus tubing has been designed specifically to fit the needs of the pharmaceutical marketplace. This silicone undergoes a battery of biocompatibility tests appropriate for pharmaceutical applications, eliminating the costly extended implant studies typically performed on biomedical grade silicones. The result: a more cost-effective silicone tubing that delivers consistent quality and the right level of assurance.

The Right Source.

Applying more than 20 years of silicone manufacturing experience, Helix Medical consistently strives to meet our number one goal: to instill confidence into every product we make. This drive to provide peak levels of assurance to our customers compels us to meet the highest standards available for our production areas.


The Right Tubing.

PharmaFocus platinum-cured silicone tubing is available in 11 popular size combinations. Typical hardness: 50 Shore A

INSIDE DIAMETER		OUTSIDE DIAMETER		WALL		ref number
inch	mm	inch	mm	inch	mm	
.125	3.18	.250	6.35	.063	1.59	61-252-21
.188	4.76	.313	7.94	.063	1.59	61-252-22
.188	4.76	.375	9.53	.094	2.38	61-252-23
.250	6.35	.438	11.11	.094	2.38	61-252-25
.250	6.35	.500	12.70	.125	3.18	61-252-26
.313	7.94	.500	12.70	.094	2.38	61-252-27
.375	9.53	.625	15.88	.125	3.18	61-252-30
.500	12.70	.750	19.05	.125	3.18	61-252-32
.625	15.88	.875	22.23	.125	3.18	61-252-33
.750	19.05	1.000	25.40	.125	3.18	61-252-35
1.000	25.40	1.500	38.10	.250	6.35	61-252-36

Values listed are typical. Please contact customer service for tolerances and custom size requests.

PHARMAFOCUS TUBING RANGE SIZE INDICATOR



Largest Size:
1.0 ID X 1.5 OD

Smallest Size:
.125 ID X .250 OD

SPECIFICATIONS

- ISO 10993-1
- ▲ FDA 21CFR 177-2600
- USP Class VI
- ▼ EP 3.1.9

50 Foot Coils (15.24 meters) • Packaged double bagged in crush-resistant boxes • Material Certification and Lot Traceability included in every box.